NitraCare Nitrile Exam Gloves

Powder-Free Non-Sterile Nitrile Examination Gloves

Characteristics and Overview:

- Finger-tip textured enable grip even when wet.
- Standard color coded sizing to identify different sizes.
- Offer excellent sensitivity for touch perception and dependable protection.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Description	Size	Packaging
NitraCare, Nitrile PF Exam Glove	X-Small	200 eaches/box 10/boxes/case
NitraCare, Nitrile PF Exam Glove	Small	200 eaches/box 10/boxes/case
NitraCare, Nitrile PF Exam Glove	Medium	200 eaches/box 10/boxes/case
NitraCare, Nitrile PF Exam Glove	Large	200 eaches/box 10/boxes/case
NitraCare, Nitrile PF Exam Glove	X-Large	200 eaches/box 10/boxes/case
NitraCare, Nitrile PF Exam Glove	2X-Large	180 eaches/box 10/boxes/case



Glove Specifications:

Brand	NitraCare
Material	Nitrile
Color	Cobalt Blue
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	XS, SM, MD, LG, XL and 2XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free, Chlorinated
Glove Exterior	Textured Finger-tips
Length (inches)	9.5 in
Palm Thickness (Avg)	3.0 mil
Finger Thickness (Avg)	3.6 mil
Tensile Strength (Mpa)	35 / 33 Mpa
Ultimate Elongation (%)	600 / 560
Pin Hole (AQL)	1.5

Quality Standards:

Manufactured in accordance with Quality System ISO 9001. Exceeds and meets current ASTM Standards: - D6319, D5151, D1671 & D6978. *Summary of Chemotherapy Testing Reports – ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifosfamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa