

# Non-Contact Infrared Forehead Thermometer

## Simple to operate

The non-contact infrared thermometer is designed for measuring body temperature, specifically for use on the forehead area. Requiring no skin contact, the thermometer works by collecting heat radiation emitting from the forehead. It's a more hygienic, reliable option for fever-checking compared to intrusive oral thermometers. Receive highly accurate measurements in seconds with the single click of a button. An integrated alarm sounds when raised body temperature is detected, making it easy to inform both user and individual of a potential fever.


- +/- 0.4° F (0.2° C)
- Auto shut-off after 30 seconds to save battery with power saver
- Pre-calibrated
- Commercial grade
- 1 year warranty
- Resolution: 0.1° F (0.1° C)




## Benefits




**FDA, CE & ISO 13485  
Certified**



**Temp Range:**  
95° F - 109.2° F  
(35° C - 42.9° C)



**Reponse Time:**  
1 second



**Testing Range:**  
2 -3 inches (5-8 cm) away  
from forehead

Patented probe design to ensure high accuracy

3 color backlight LCD for easy reading

Storage of the last 10 readings

°C and °F readings available



#### Forehead/Human Mode

Light Indicator	Temp Range
Green	95° F - 99.3° F
Orange	99.5° F - 100.2° F
Red	100.4° F - 109.2° F ( <i>"HI" will appear if exceeded</i> )

Prod No	Description	Color	Reponse Time	Accuracy	Temperature Range	Pack
THDG986	Non-contact IR Forehead Thermometer	White	1 second	+/- 0.4°F (0.2° C)	95° F - 109.2° F (35° C - 42.9°C)	1 each (60 carton)

**San Jamar**  
555 Koopman Lane  
Elkhorn, Wisconsin 53121  
USA  
T: +1.262.723.6133  
F: +1.262.226.8664  
info@sanjamar.com  
www.sanjamar.com

**Canada**  
3280 Bloor Street West  
Suite 1140  
11th Floor, Centre Tower  
Toronto, ON M8X 2X3  
CANADA  
T: +1.262.723.6133  
F: +1.262.226.8664  
canada@sanjamar.com

**Europe**  
Schoorstraat 26a, bus 1  
2220 Heist-op-den-Berg  
BELGIUM  
T: +32 15 22 81 40  
F: +32 15 22 81 48  
emea@sanjamar.com

**México**  
Av. Paseo de la Reforma No. 350 - 10  
Col. Juárez  
Deleg. Cuauhtemoc, CP 06600  
MEXICO, D.F.  
T: +52 (55) 3626 0772  
mexico@sanjamar.com

 **san jamar**<sup>®</sup>  
smart. safe. sanitary.